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English

These instructions apply to all the versions of the OTOPLUS unit and all the related accessories. For this reason, not all the paragraphs of this manual will refer to the equipment you have bought.

MEDI-CARE SOLUTIONS will provide the authorised service centre with all the information necessary to install and repair the equipment.

No part of this manual is to be reproduced, stored in a retrieval system or transmitted in any form or in any means, i.e. electronic, mechanical, photocopying, translation or otherwise, without the prior written permission of MEDI-CARE SOLUTIONS.

The OTOPLUS unit is a complex for otolaryngology whose use is permitted exclusively to doctors or specialized in ENT field.

The device in object must be used matching (electrical connection) just other medical devices **and with exclusive use destination the patient area**

The following MD should not be used near or on top of other equipment, if necessary, use it near or overlap with other units to observe the DM and control normal operation in the configuration that is commonly used.

The expected useful life of the device is 10 years

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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- 1.1 Intended and environment of use
 - 1.1.1 General characteristics and identification of the parts
- 1.2 Cautions
- 1.3 Symbols

2. INSTALLATION**3. USE OF THE INSTRUMENTS**

- 3.1 Front commands
- 3.2 Insufflator
- 3.3 Syringe
- 3.4 Suction
 - 3.4.1 Suction accessories
 - 3.4.2 Antibacterial filter
- 3.5 Medicine nebulizer

4. USE OF THE ACCESSORIES

- 4.1 Pole for microscope and monitor support
- 4.2 Germicide compartment
- 4.3 Clar lamp
- 4.4 Larynx mirror-warmer
- 4.5 Optical housing support

5. MAINTENANCE

- 5.1 Overall control
- 5.2 Filling the water tank
- 5.3 General maintenance and spare parts
- 5.4 Cleaning of the unit and secretions tank
- 5.5 Single use secretion tank
- 5.6 Empty condense filter
- 5.7 Discharge multipurpose secretion tank
- 5.8 Disinfectant products chart

6. TECHNICAL CHARACTERISTICS

- 6.1 Technical data charts – absorbed power
- 6.2 Working periods of the instrument summary – identification of the unit

7. IDENTIFICATION LABELS


- 7.1 Label data
- 7.2 Fuses identification label
 - 7.2.1 Cables and sectionig tube
- 7.3 Disposal rules

8. GUIDE CHARTS AND PRODUCER'S DECLARATIO

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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1. INTRODUCTION

1.1 Intended of use and environment of use

The devices in question are intended to carry out treatments and/or specialist outpatient visits in the ENT sector or affecting the ear, nose, throat with the possibility of carrying out cleaning and/or inspection operations on the patient's auricle using active functional elements and /or to carry out targeted procedures to verify physiological functions in the ENT field. The operators assigned to use the devices must be specialized (medical professionals specialized in ENT) or duly trained (clinical nurses or operators specialized in the ENT sector). It is not envisaged that the dental units can be used for activities other than the ENT sector and by non-specialised (lay people) or duly trained personnel.

The intended environment of use is exclusively in the ENT specialist medical sector (hospitals, clinics, polyclinics, private and public specialist medical practices in the ENT sector, etc.), compliant with the electrical directives in force and specific for medical environments; It is not intended for use in different environments (e.g.: for home use, non-medical environments, veterinary environments, mobile units (ambulances, etc.)).

1.1.1 General characteristic and identification of the parts

Unit for otorhinolaryngology completely autonomous and mobile on castor wheels composed of (see pic. 1):

- Instrument tray holder, which is able to hold up to a maximum of 5 instruments. It has articulated arm which enable one to position it in the most conformable way for the professional's activity. On the board is placed the electronic console with LCD display and touch screen orientable.
- Modular cabinet containing the functional equipment of the unit
- Anti-scratch working top with drop-saving back border embodied in the top; or alternatively desk for instruments trays with plastic cover.
- Housing for the auxiliary instruments complete with a series of drawers & spaces of various sizes and dimensions according to one's choice and microscope and monitor pole support.

1.2 CAUTIONS

- To move the unit with the drawers closed and pc-monitor and Microscope within rest position. Unit weight it's about 180kg (based on configuration).
- Do not install the machinery in environments which are exposed to explosions or frost
- The temperature of the environment in which the equipment is installed and used must be between 10 °C and 30 °C. The relative humidity environment of use or storage must be understood between 30% and 70%; pressure 500hpa at 1060hpa
- During installing, check that the electrical tension on the sticker on the machine corresponds to that in use in order to avoid damaging.
- At the time of installation, the correspondence between what is indicated on the plate on the back of the equipment and the voltage in use must be verified, in order to avoid damage.
- The building, repairing, modifying, calibration and all the other operations that require personnel authorized by MEDI-CARE SOLUTIONS should carry out the opening of the unit.
- The electric circuit in which the unit is connected to, should comply to the IEC 60364-7-710 regulation (regulation regarding electrical circuits in rooms used for medical purposes).
- The equipment should be used exclusively following the instructions illustrated in this manual.
- 9) MEDI-CARE SOLUTIONS refuses any responsibility for the following: damages or accidents caused by the carelessness of the user, the non-compliance to the safety rules, a non-correct use of the machine, a use of the machine which is different to that of its purpose. MEDI-CARE SOLUTIONS denies all responsibility which has nothing to do with faults in the material and non-correct building.

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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10) Before installation and functioning check that there are not loose parts, burnt electric cables, or various anomalies and if they are evident, please contact Medi-Care Solutions or authorised technical staff.

11) Mobile phones and radios can influence the good functioning of medical device OTOPLUS.

12) If during the use of unit OTOPLUS happen strange changes of performance or various anomalies we suggest to interrupt all operations, switch off the machine and contact authorised staff Medi-Care Solutions in order to avoid any kind of damage.

13) Do not install this unit in places exposed to the risk of explosions, rich in oxygen or exposed to frost

14) The DM is not protected against splashing water and foreign bodies: **DEGREE IP20**

15) **WARNING!** To avoid the risk of electric shock, this equipment must only be connected to a supply main with protective earth."

16) **WARNING!** There shall be no modification of this device

17) **WARNING!** The MD in question should not be used adjacent to or stacked with other equipment, if necessary, use it close to or stacked with other equipment it is necessary to observe the DM to check whether, in the configuration in which it is used, its operation remains normal

18) The MD can be transported, stored with the environment:

- Temperature: da -20°C a 60°C;
- Humidity: 10% - 90% ;
- Pressure = by 500hpa at 1060hpa

19) The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment"

20) **IMPORTANT!** In the event of a serious accident occurring in relation to the device, it is necessary to report to the Manufacturer and the competent authority of your country.

• **APPLIED PART OF THE UNIT**

1. Terminal of Auricular wash syringe
2. Suction terminal (cannula and funnel)
3. Insufflation Terminal (cannula and disposable terminal)

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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1.3 Symbols

Classification of the "BF" type machinery. "BF" type machinery (machinery of the type I,II, III or machinery housing an internal electric source with a protection degree against direct and indirect contact and in particular for:

- permitted current dispersion
- reliability of the ground contact)

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This symbol indicates: TO CLARIFICATION CONSULT THE ATTACHED DOCUMENTATION

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This symbol indicates: ATTENTION!!! CONSULT THE ATTACHED DOCUMENTATION BEFORE USE THE UNIT

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Symbol that indicates "generic warning".

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Earthling contacts of the machinery

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**1.4 WARNING, POSSIBLE SIDE EFFECTS**

Always remember, before you use the MD (medical device) to make a careful medical history and a careful examination of the patient to recognize situations in which it recommended the use of active instruments such as the ear wash syringe, suction, insufflation, that can cause side effects such as redness, bleeding, abrasions, minor burns, etc..

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**CONTRAINDICATIONS**

The contraindications are those specific situations in which the use of the tools are preferably not active carried and consequently the device must be used where appropriate with a restriction defined relative to the cost / benefit ratio resulting from the use of the system itself. It is therefore considered that any adverse or residual risk can be defined by MEDI-CARE; is the responsibility of the specialist do the analysis costs / benefits of the use of MD based on their knowledge and experiences.

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**2. INSTALLATION**







The installation of the equipment must be made respecting the instructions of the following booklet of instructions. It must be performed by authorized personnel MEDI-CARE SOLUTIONS or, given the simplicity of the operations can also be performed by the customer under its responsibility.

The user should always make sure that the electric circuit complies to IEC 60364-7-710 regulation or special rules of your country.

The operations are (after unpacking the unit and take IFU):

- Taking unit out from the pallet from basement or solid parts and not from working top or instrument tablet,
- Unpacking the instruments and put them on the tablet as picture 1
- Unblock the wheels brake and positioning unit
- Screw, if it's present, the monitor and microscope support pole
- Full the tank but not to the rim (See. Chapter 5.1)
- Insert the plug in a socket ➔ It is advisable to insert the plug in a socket with a bipolar switch compliant to IEC and having the following technical characteristics:
 - 250 Vac - 50 A contact opening of 3 mm

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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3. USE OF THE INSTRUMENTS

In order to use the instruments of the machine it is necessary to proceed with the turning on of the machine by pressing the switch on front of the work bench (see pic.1). When the green light is on, this means that the machine is working. Turn it off by pressing the same switch again.

3.1 Frontal commands

The following commands are present on the front panel of the unit: it is possible to use them only after turning the general switch on (see picture 1)

- 1- General switch
- 2- Water heating switch
- 3- Potentiometer for clar lamp. By turning it clockwise the light increases by turning it the opposite way it decreases

3.2 Insufflator (Picture 2/a) (optional)

For a correct use of the instrument, it is advisable to work according to the following periods: Work 5 min - Off 10 min. The instrument delivers compressed air up to a maximum of 2.2 Bar for insufflation within ENT.

The functioning of the instrument starts when it is taken out of its housing and pushing the black pushbutton. When the corresponding light on the instrument holder is on, this means that the instrument is working. In order to set the intensity of air, one should adjust it with the tap under the tablet.

WARNING! In the event that the activation of the instrument, the compressor does not work, press the black button to download the air in the circuit; see Pic.3.2

Notice: All these parts are washable or chemically sterilized every time you want according to their internal procedures. See point 5.8 table B.

Notice: The cannula is a medical device; for cleaning and sterilization procedure see user manual of manufacturer.

Notice: The supplied nozzle is connected on the cannula; it is disposable therefore it must be replaced after each use/application. Notice: Use only original spare parts.

Fig 3.1/Pic. 3.1

3.3 Syringe for auricular washing (Picture 2/b) (optional)

For a correct functioning of the instrument, we suggest to use it according to the following periods:

Work 5 min → off 10 min

The instrument provides approximately 350 ml / min of deionized water 37°C (±1°C) at a maximum pressure of 2,2 bar. Before using the syringe, it is necessary to fill up, not completely full, the tank.

The functioning of the instrument is activated by taking it out of its housing. When the corresponding indicator-light on the instruments tablet is lit up, this means that the instrument is functioning. For obtaining warm water flux, firstly turn on the warm water switch placed in the instrument tablet of the unit (see section 3.1): the led stand-by on the tablet starts blinking. After some minutes when the right temperature is reached, the led is lit up and the acoustic sign is activated. Now it is possible to push one of the two buttons placed in the syringe hand-piece.

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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! NOTICE: The first heating cycle can last some minutes (about 10 minutes).

! Attention! In the event of a reserve water an intermittent audible alarm appears. Turn off the machine and refill the tank as 5.1 chapter.

! Notice: After a prolonged pause with the machine off, when the temperature is reached, it is always advisable to let come out from the syringe a water flow for about 10 seconds so that to asset the machine.

! Warning! In the evidence that the water temperature overcomes the set value, a safety thermostat blocks automatically the washings syringe functioning, until the normal cycle is re-activated. In that case appears the intermittent sound to indicate that the thermostat is active. In that case, it is necessary to contact an authorised technician of MEDI-CARE SOLUTIONS, which is the only authorised to intervene.

! Note: The syringe is a medical device certified by Luzzani dental; for his maintenance and sterilization, please follow the instructions in the Technical Manual here attached

! Notice: Use only spare parts Medi-Care Solutions.

3.4 Suction (See Pic. 2/c)

For a correct functioning of the instruments, we suggest to respect the following using times: **Work: 5' - Pause: 10'**

The instrument allows you to suction with a maximum vacuum till -0.85 bar $\pm 5\%$.

The instrument's functioning is obtained by extracting it from its site. On the tablet a green LED will be switched on once the instrument is extracted. Turning the knob to the right the suction power is increased and vice versa).

! Notice: If the container is full of secret operative system overflow aspiration is stopped, so to drain (see 5.2)

! Notice: The useful capacity is stated on the vessel, the meeting cannot work on inclined planes.

! Notice: In case of dry pump, in order to improve performance, after you put the suction in the its support, the pump remains active for about a minute.

! NB: Use only spare parts by Medi-Care Solutions.

Fig 3.1a/Pic. 3.1a

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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**TUBO PROVENIENTE DAL RIUNITO
TUBE COMING FROM THE UNIT**

3.4.1 SUCTION ACCESSORIES

The unit is set with a pulling-out neck headed by a piece of silicone pipe (see pic. 3.4.1) on which it is possible to connect the suction iron cannulas; moreover, together with the machine is supplied a kit for setting the funnel for secretions collecting (see pic. 3.4.2) which contains a pipe fitting to be coupled to the suction pipe + funnel (see assembling phase).

Notice: All these parts can be cleaned or washable according to their internal procedures with common disinfectants (See table 5.8); except suction tube the parts they are autoclavable (See table 5.8).

Fig 3.4.2
Pic. 3.4.2

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NOTICE: The capacity of the secretion tank are applicable to 1 and 2 liters.

WARNING: Do not keep the rolling shutter always closed, because it could cause suction locking.

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3.4.2 Antimicrobial hydrophobic filter

On the Unit inside the side door is planned an antimicrobial hydrophobic filter on suction, in order to ensure the killing (99%) of any bacteria that come from the suction system. The replacement **is advisable** every 2-5 working days.

WARNING! To change the filter, follow the direction indicated on the filter (IN) and pipes.

Notice: Medi-Care Solutions recommend the use of pvc pipes or frosted silicone size of 8x12mm

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Supporto ampolle/Ampoule holder

**Fig. 3.7
Pic. 3.7**

3.5. Air syringe for atomizer and insufflator (opt.) (Pic. 2/d)

For a correct functioning of the instruments, we suggest to respect the following using times: Work: 5 min / Pause: 10 min.

This instrument provides dry compressed air to a maximum pressure of 1.8-2 bar; at the gun you can connect the atomizers or a choice of air blow through the politzer cannula. The functioning of the instrument is obtained by extracting it from its seat (see pic. 3.7). On the instruments tablet will light on a green LED that advising you are using the instrument. After the application of the ampoule to the syringe to obtain the vaporization of the medicaments products, operate on the lever of the instrument hand piece (see pic. 3.5).

It is possible to adjust the pressure and the flow of air through the tap under the tablet. In counter clockwise to increase the pressure and vice versa (see pic 3.6)

WARNING! To avoid keeping the circuit too much under pressure, we suggest not to keep the syringe extracted too long without using it.

NOTICE: USE ONLY SPARE PARTS supplied by MEDI-CARE SOLUTIONS.

Warning: For maintenance, cleaning and sterilization, follow the table 5.8 attached to this document

**Fig. 3.6
Pic. 3.6**

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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4. USE OF THE ACCESSORIES

4.1 Pole for microscope and monitor support (pic. 4/a)

Optionally, a pole can be arranged on the top (to support a monitor or tray holder or a microscope). In this case the pole is already mounted on the unit type; see photo 4.1.

Attention! Installation of the monitor must be carried out by specialized personnel and by participating in the manufacturer's user manual.

Attention! It is recommended not to move the unit when pc-monitor is rotated backwards the unit.

4.2 Germicide compartments (See Pic.4/b)

On request, the machine can be equipped with this auxiliary drawer that permits only the maintenance in hygienically conditions of the instruments already sterilized.

The functioning of the UV lamp is obtained by turning the general switch of the unit on (see chapter 3).

The germicide lamp can be turned on and off by a switch placed on the side of the drawers turns the lamp off automatically if the germicide compartment is accidentally opened, and it will turn the lamp on again when the drawer is closed.

Attention! We repeat once more that the germicidal space does not constitute substitutive element of the sterilization process, but function only as co-adjutant and supporter.

Notice: Do not use lamps with technical features different from those indicated by Medi-Care Solutions.

Notice: We suggest to substitute the lamp each 500 hours of functioning

ATTENTION! Take the plug of the Unit out before carrying out any operations on UV circuit.

4.3 Clar Lamp (Pic. 4/c)

The lamp is turned on automatically by taking it from its housing, to turn it off it is necessary to place it back.

The adjusting of the intensity of the light is obtained by adjusting the potentiometer placed on the front (see chapter 3.1)

4.4 Mirror-warmed larynx (Pic. 4/c)

The functioning of the mirror-warmed larynx is obtained by pressing the switch placed next to the fan (see figure).

Place the surgical instrument in front of the warm air flux and wait a few seconds.

WARNING! When overused, the system has a thermal safety switch which stops the air flow until normal values of temperature are reached again; it is necessary in that case to wait for some minutes before using it again.

4.5 Cleaned and Dirty Optics Housing (pic. 4/d)

On the Unit is standard planned an integrated support for the dirty and cleaned optics. There are 2 closed supports for inserting some disinfectant solutions and for the resting optics. The supports is longer 250mm for rigid endoscope and 500mm for flexible endoscope.



Warning! Do not use aggressive solutions and do not maintain them on the support more than 24 hours; wash the supports with hot water after use.

Warning! To preserve the supports we suggest using solutions for the sterilization at cold or generators of free in order to conserve the supports

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

Vano Germicida/Germicide Compartments
Accensione Cassetto UV/UV lamp power on
Accensione Scaldaspecchiotti/ Power on mirrors heating system

Gancio per supporto e accensione lampada frontale Support Hook and power on Head Lamp
SUPPORTO OTTICHE
Fixed Holder for endoscope

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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5. MAINTENANCE

Warning! The following proceeding have to be carried out when the machine is not plugged into the mains and with specific PPE aids to avoid biological contamination.

5.1 Empty of the condense filter

It is advisable to check the conditions of the containers and of the accessible internal parts cyclically

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5.2 Filling the water tank

In acoustic case of signal intermittent it means that it is intervened the reserve of the tank and that it is necessary to fill it up. In order to fill the tan please proceed in the following way:

- 1- Take the plug out of the machine
- 2- Open the door of the auxiliary instruments compartment
- 3- Take the tank out, unwind the tap and with a funnel pour the water, without reach the maximum level.
- 4- Close the tap and the door
- 5- Plug the machine in again


Notice: It is recommended to use only distilled H2O mixed with disinfectants solutions to avoid problems of limestone deposits in the tubes and proliferation of bacterial charges.

WARNING! Do not leave water in the tank for more than 12 hours; in this case, to avoid bacterial contamination, it is recommended to empty the tank and dry it with a sterile cloth or alternatively, add a commercially available preservative solution to the water

Aprire lo sportello/
Open the door

Vano alloggiamento
Auxiliary instruments
compartment

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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
Rimuovere pannello/To reach the pump take off the rear panel.

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5.3 General maintenance and spare parts


To guarantee the correct functioning of the instrument Medi-Care Solutions suggests to effectuate a periodical revision of the entire machine.

Medi-Care Solutions does not take charge of any responsibility for damages, hitches or inefficiencies due to the substitution of non-original spare parts.

5.3.1 OIL PUMP MAINTENANCE (optional)

1. Oil Level Check

It is suggested to check the oil level of the pump's engine every 2-3 working months, according to the hours of activity sustained by the Pump (See pic. 1)

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Fig. 1
Pic. 1

2. Oil Refill

After approximately 500h of suction use, and in case of darkoil due to its deterioration, proceed to its replacement by following point 5.4 of instructions manual DVP (see pic. 1).

3. Filter Replacement

The pump's filter is able to reduce residual fumes for 2000 working hours, after such period of time it is necessary to replace it with a new filter; see the following pictures 2, 3:


Fig. 3
Pic. 3

B. svitare il filtro/To unscrew the old filter and to replace it with a new one.
Attention: Non forzare il filtro/To tighten gradually without excessive force.

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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**5.4 Cleaning the unit**

The choice of materials used for building the unit has been carefully taken in order to correspond to the cleaning and disinfecting needs. Given the big differences between chemical products used in a medical study, it is possible to ruin the tops. The major damages are caused by the amount of time that these chemical products are left on the tops. It is important to immediately dry the interested parts with a damp cloth.

To clean the painted or metal parts, use a specific product for these materials which doesn't damage their composition.

IMPORTANT! Do not use products with alcohol, ammonia or petrol

5.4.1 Cleaning and maintenance of secretion tank

The tank, which contains the aspirated liquids, can be cleaned with water and detergents, which are compatible with the materials (see chart) used and/or sterilized in an autoclave as of manufacturer instruction of use.

Separate the main parts of the tank before the sterilization insert the glass upside down in the autoclave and do not place anything on top of it whilst disinfecting.

DO NOT USE SOLVENTS OR ALCOHOL FOR CLEANING AND DISINFECTING. THE USE OF THESE PRODUCTS CAN DAMAGE THE PLASTIC PIECES OF THE SYSTEM.

Always substitute the damaged or faulty parts with original spare parts supplied by MEDI-CARE SOLUTIONS.

Notice: MAK's suction jar using by Medi-Care are certified medical devices; for maintenance and sterilization, follow the manufacturer's instructions attached to the following.

5.5 Disposable secretion jar

The unit can be set for single-use kit of different type. **Please keep attending during the replacement of the ampoules in order to avoid inverting the suction tubes VACUUM and PATIENT, because it could damage the motor.**

Notice: If the intake of liquids or solids into the engine, immediately turn off the suction and call the service center for maintenance.

Notice: Serres's suction jar using by Medi-Care are CE certified medical devices; for maintenance and sterilization, follow the manufacturer's instructions attached to the following

5.6 Empty of the condense

Near the panel there is a filter for collecting the condense created by the compressor; in case it is filled it is necessary to empty it; place a glass under the filter, turn the knob and pull downside the knob, empty the condense and close again the knob!

Warning! If the knob is not closed the compressor leaks air through it!

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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5.7 Empty of multipurpose secretion tank

At the end of a day's work or when the overflow system (valve included in the jar) is running, so suction is forbidden it is necessary to clean the tank where all the secretions gather. In order to carry out this operation, proceed in the following way:

- 1- Unplug the unit
- 2- Open the door where the tank is housed
- 3- Unwind the tank by turning it anticlockwise to separate it from the cover which is fixed
- 4- Empty the tank and clean it thoroughly
- 5- Re-screw the tank to the cover, close the door and plug the unit again.

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WARNING! It's not recommended to aspirate large amounts of fluids or disinfectant to prevent it from being drawn into the circuit by stopping the pump; we suggest engaging the suction cannula with a maximum diameter of 2/3 mm and aspirate the disinfectant sure it does not may trigger the gurgling or splashing in the secretion tank.

Notice: MAK's suction jar using by Medi-Care are certified medical devices; for maintenance and sterilization, follow the manufacturer's instructions attached to the following

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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5.8 Cleaning and disinfection/sterilization of the part

A. Instrument and detachable parts sterilization chart (if present)

INSTRUMENT/PART	SUPPLIERS	KIND OF STERILIZATION
Ear Wash Syringe (terminal, body)	LUZZANI	See user manual of manufacturer
Iron suction cannula	GIMA	See user manual of manufacturer
Plastic cannula	CATTANI	See user manual of manufacturer
Silicon Suction Pipe 6x12	GIMA	See user manual of manufacturer
Canister multipurpose (MAK series)	FLOWMETER	See user manual of manufacturer
Disposable bag container (Serres)	SERRES	See user manual of manufacturer

B. Cleaning and disinfecting products chart external parts

WARNING! Do not use different products without first contacting Medi-Care Solutions

Disinfectants	Pharmaceutic products and cosmetics
+ bactol 5% + clorammina - DDT + delegol 5% = dimammina 5% = iodine dye - fenic acid + lysoform 2% - TB lysoform + bactol + merfen 2% + octozone + perydrol + resorcina 1% = sagrotane 5% + alcohol to be purely burnt + mercury chloride + trosilina G + oxiginated water = zefirol + sodium ipocloride (amuchina) + free radicals (ec-ster)	+ blood plasma + sun tan lotion + hydroplex = iodine dye + lanolyn = mantolo 90% in alcohol - nail varnish - solvent + mouth wash + periston R + vaselyn + vicks rub
CHART SIGNS + compatible - non compatible = compatible but with caution	

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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6. TECHNICAL FEATURES

6.1 Technical data chart

Electrical tension	230VAC +/-10%
Frequency	50/60 Hz
Power absorbed by the unit (according to the configuration)	1200 VA
Weight	120 Kg
Maximum number of instruments	4 pcs
Maximum air pressure	2.8 Atm
Air capacity	26 l/m
Maximum water pressure	2.2 Atm
Water capacity	0.4 l/m
Dimensions (in mm):	
Height of the working bench	805 mm
Width	860 mm
Depth	745 mm

Absorbed power chart

DESCRIPT.	ABSORB.	TENSION	SERVICE
Germicide Lamp	11 W	230 V ac	Continuos
Aspirator	120 W	230 V ac	Continuos
Compressor	120 W	230 V ac	Continuos
Pump	12 W	230 V ac	Continuos
Electrovalve	24 W	24 V dc	Continuos
Resistance	800W	230 V ac	Discontinuos

6.2 Summary of the working periods of the instruments

INSTRUMENT	WORK	OFF
Syringe	5'	10'
Insufflator	5'	10'
Suction	5'	10'
Atomizer	5'	10'

1	2	3	4	5	6	23	24
7	8	9	10	11	12	26	25
N° matricola							

Building Datas

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

6.3 SCHEDA TECNICA/TECHNICAL FEATURES

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NOTE: MEDI-CARE SOLUTIONS on request will provide electrical circuit diagrams, the parts list, descriptions, calibration instructions or other information to assist the AUTHORIZED SERVICE PERSONNEL in the repair of parts of the MD that are considered repairable.

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The Serial number **label** is located in the back part of the unit, and the details are follows:

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In the back part of the unit is located the label of the fuses – its values are listed in the chart 6.2 of the following MI.

Inside the machine, the protective mass is indicated by the following stucked symbol: 

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NOTE: THE FUSES OF THE DEVICE ARE IN THE BACKSIDE TAKING OFF THE FIRST DRAWER.

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www.euroclinic.it; info@euroclinic.it.

7.3 DISPOSAL

According to the legislative decree 14/03/2014, n. 49, "Implementation of the European Directive 2012/19/UE about decrease of dangerous substances on electric and electronic equipment" and moreover it refers to garbage disposal, the symbol of crossed out garbage bin (on the back of the equipment) indicates that the product at the end of its working life must be collected separately from other garbage according to the following law 31 may 2016 n° 121 and 8 march 2010, n.65 or follow your country laws.

The separate collection of rubbish of this equipment (in disuse) must be arranged and coordinated by the manufacturer. When the user wants to throw away the equipment, he has to contact the Manufacturer Company that will indicate all procedures to follow for the disused equipment separate collection of rubbish.

Therefore, the appropriate separate collection of rubbish, aimed to actuate afterwards disused equipment recycling, treatment and to ambient garbage disposal, contributes to avoid possible bad effects on the environment and on the health. In addition, these procedures raise the reemployment and/or the recycling of the materials of which the equipment is composed.

The unauthorized garbage disposal, on the part of the holder, requires administrative law sanctions application.

MEDI-CARE SOLUTIONS S.R.L. remains available for further information about this matter:

Registered Office: MEDI-CARE SOLUTIONS Ltd

Via Della Zecca 1, 40121 Bologna (BO)

Facility: Via Pietro Nenni 3, 40026 Imola (BO) - ITALY

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The electro-medical device requires special precautions regarding EMC (Electromagnetic Compatibility) and must be installed in accordance with the EMC information contained in the accompanying documents.

Portable and mobile radio-communication equipment can affect the operation of EUT (Equipment under test).

The use of accessories, transducers and cables other than those specified, with the exception of those sold by the manufacturer as spare parts for internal components, may cause an increase in emissions and a decrease in immunity.

The device or the system must not be used near other devices and, if it is necessary to use it near other appliances, the electro-medical device must be observed to check normal functioning in the configuration in which it is used.

The electro-medical device has been tested and founded in compliance with the emission and immunity limits of electromedical equipment according to IEC60601-1-2: 2014. These limits are specifically designed to provide adequate protection against dangerous interferences in a typical medical room. However, there is no guarantee that interferences will not occur in a particular installation. If the electro-medical device, interacting with another one, causes or receives detectable interferences, the user is invited to limit them by adopting one or more of the following measures:

1. Reorient or relocate the receiving device;
2. Increase the distance between the appliances;
3. Connect the equipment to an outlet on a circuit different from the device (s) causing the interference;
4. Contact the manufacturer or local technician for assistance.

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data

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Proximity magnetic fields

IEC 61000-4-39

Test frequency	Modulation	IMMUNITY TEST LEVEL (A/m)
30 kHz ^{a)}	CW	8
134,2 kHz	Pulse modulation ^{b)} 2,1 kHz	65 ^{c)}
13,56 MHz	Pulse modulation ^{b)} 50 kHz	7,5 ^{c)}
a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) r.m.s., before modulation is applied.		

12	Francesco Marrone	11-11-2023	DG Fabio Casella	11-11-2023
Release	Redatto	Data	Approvato	Data